


Relipore Xtreme

Advanced Sterile Adherent Dressing Pads



Why are these dressings unique?


- Advanced wound contact layer.
- Low-adherence technology.
- Hi-tech hexagonal mesh.
- Fluid flows in one direction only.
- Absorbs up to 10 times their own weight.
- Retains fibre integrity when cut.
- For use on lacerations, minor burns, grazes and as a post-operative dressing.



| REF | PRODUCT | SIZE | UNIT | PACK |
|------|--|------------|-----------|------|
| 2600 | Relipore Xtreme Adhesive Dressing Pads Sterile | 8cm x 10cm | Box of 50 | 1 |
| 2601 | Relipore Xtreme Adhesive Dressing Pads Sterile | 8cm x 16cm | Box of 50 | 1 |




ReliporeXtreme - advanced adherent dressing pad.



Hi-tech hexagonal mesh allows fluid to flow in one direction only.

Pad absorbs up to 10 times its own weight and retains its fibre integrity when cut to size.

Flexible, breathable stretch fabric conforms to the body.



The fibres in the absorbent layer are arranged horizontally so that they encourage fluid to spread evenly across the full breadth of the pad.

ADVANCED WOUND CONTACT LAYER

ReliporeXtreme uses an advanced wound contact layer that encourages fluid to pass from the surface through to the absorbent layer, but helps prevent it striking back through. This is achieved using hexagonal shaped pores that have a burr on the underside of the layer. The layer is constructed from a nylon polymer that is designed not to stick to a wound site.

The ReliporeXtreme dressing has the wound contact layer on one side - and on the other - a flexible, breathable, adherent stretch fabric that conforms to the body, is easy to remove, but remains securely in place once applied.

HIGHLY ABSORBENT LAYER

ReliporeXtreme has 350g /m² of advanced absorbent material, over 3 times the weight of standard Relipore! This absorbent material is capable of absorbing up to 10 times its own weight. The fibres in the absorbent layer are arranged horizontally so that they encourage fluid to spread evenly across the full breadth of the pad, maximising its absorbent capability and preventing strike through and strike back.

PACKAGING

The pads are packaged in an easy to open pouch that can act as a sterile field for the dressing whilst a wound site is being prepared. One side of the pouch is transparent, assisting clear product identification in low light conditions.

CUTTING RELIPORE

Relipore retains its fibre integrity when cut to size, thus maximising its flexibility.

INDICATION

ReliporeXtreme is indicated for use on lacerations, minor burns, grazes and as a post-operative dressing. It is equally effective on low, medium and high-level fluid discharge situations.

CE CLASSIFICATION

Relipore is a Class 1 Sterile Medical Device.

Advanced wound contact layer.

Low-adherence technology.

Hi-tech hexagonal mesh.

Fluid flows in one direction only.

Absorbs up to 10 times their own weight.

Retains fibre integrity when cut.

For use on lacerations, minor burns, grazes and as a post-operative dressing.

