

TraumaFix® Dressing

Designed for major bleed trauma

PATENT APPLIED FOR

IMPROVED NOW WITH PRESSURE PILLOW



Innovative new pillow increases the pressure to stem major blood flow. Major bleed technology. Amazing fluid absorption.

What do you have in your First Aid Kit to deal with a major bleed injury?

Ambulance dressings have their place, but they are not substantial enough to cope with a major trauma bleed, and applying enough pressure with a light conforming bandage will be near impossible. Abdo-pads might be substantial enough to cope with the bleed, but they need a suitable pressure bandage to attach them, which then needs to be secured. This is very time consuming, difficult with one injury and dangerous with multiple bleeds.

All too often, non-sterile sanitary towels are used as major bleed dressing pads. Typically these contain super-absorbent powders that clot blood into a gel on contact. Just one grain of this powder has the capacity to cause a deadly blood clot if it were to ingress into the casualty's venous system.

TraumaFix® offers a substantial, highly absorbent dressing pad, capable of holding over ten times its own weight in fluid. A strong bandage allows direct pressure to be applied quickly and maintained easily. This is secured instantly with a Velcro® style fastening.

Multiple major bleed injuries can be stemmed quickly, allowing the casualty to be moved easily. The dressing is a significant advance in the treatment of trauma injuries.



Place highly absorbent dressing pad with pressure pillow over trauma...



Direct pressure is applied with 4.5m of pressure bandage...



The bandage is secured instantly with a velcro style fastening...



Major bleeding is stemmed quickly, allowing the casualty to be moved easily.

REF	PRODUCT	DESCRIPTION	SIZE	UNIT	MOQ
962	TraumaFix Dressing	Sterile	10cm x 18cm	1	10
963	TraumaFix Dressing	Sterile	15cm x 18cm	1	10

