

### Data Sheet For The Cleansing Wipe

No.	Item	Unit	Standards	
<b>Package bags pouch</b>				
1	Appearance	Surface Quality	/	The Surface Should Be Uniform, No Wrinkles, Cuts, Sand Holes, Holes, Hard Blocks Or Paper Layer Adhesion. Appearance of the printing clear, no ghosting
		Laminating Quality	/	Inner membrane should be smooth, flat, uniform and solid
2	Tensile Strength (CD)		KN/m	≥2.000
3	Peeling Strength	Paper&PE	/	Un strippable
		PE&Aluminum Foil	/	Un strippable
<b>Non-woven pad</b>				
4	Gram Square Meter		g/m <sup>3</sup>	49.8
5	Thickness		mm	0.34
6	Size		mm	100(±2)*80(±2)
7	Breaking Strength	CD	N/5cm	30.2
		MD	N/5cm	86.5
8	Elongation at break	CD	%	146
		MD	%	38.6
9	Fluorescent material		/	No
<b>Liquid</b>				
10	Concentration Of Sodium Chloride		%	0.85-0.95
11	Liquid content		ml	2.0(±0.1)