


<b>Product Name</b>	<b>Non woven tape</b>	
		
Specification	616 617	1.25cmx10m 2.5cmx10m
Technical Data	Material	white non woven 22g/sqm imported wood pulp and PVA as raw material, made by a special process to wet nonwoven fabrics
	Adhesive	hot-melt adhesive/ Henkel Synthetic rubbers           20 – 30% Hydrocarbon resins       50 – 70% White mineral oils         10 – 20% Antioxidants                 0.1 – 1.0%
	Coating Weight	30±2g