

Beat CPR™



CPR Resuscitation Aid



CPR Resuscitation Aid

The Beat CPR device provides first aiders with a clear, audible tone, helping maintain the compression rhythm at optimum efficiency. The device is pre-set at 110 beats per minute, which may be increased or reduced by the user.

1. LCD Indicator
2. Power
3. Parameter Adjustment
4. LED Rhythm Indicator
5. Speaker



- Compact.
- Audible tone guide CPR compression speed.
- Pre-programmed to 110 beats per minute, adjustable.
- Display visible in all light conditions.
- Long-life renewable battery included.

REF	PRODUCT	UNIT	MOQ
2860	Beat CPR™ - CPR Resuscitation Aid	1	1

