

Advanced Fabric Range

Now lighter-weight, made from modern hi-tech and extra flexible fabrics that are more comfortable, conformable, washproof and easy to wear.

The elastic fabric stretches in all directions, ideal for covering joints and skin folds with ease and comfort. The fabric coated with a hydrophilic polymer making it water repellent and washproof. The new super-sticky zinc oxide adhesive ensures that Dependaplast Fabric Plasters remain stuck firmly in place. The adhesive comes cleanly away, avoiding unsightly plaster marks being left on the skin.

Dependaplast plasters offer a full and comprehensive range of the most popular shapes and sizes. Plasters are now neatly packed in a variety of boxes, allowing for clear, easy access, and because the boxes are smaller, save on storage too. Ideal for refilling first aid kits or to keep as backup supplies. All Dependaplast plasters are latex-free with a Rayon wound contact layer.



REF	SIZE	BOX QTY	MOQ
555	2cm x 12cm	50	1
509	4cm x 2cm	100	1
511	4cm x 4cm	100	1
512	6cm x 2cm	100	1
514	7.5cm x 2.5cm	100	1
515	7.5cm x 5cm	50	1
513	7cm x 2cm	100	1
519	Anchor	50	1
518	Fingertip	50	1
516	Assorted 30 4cm x 4cm 30 7cm x 2cm 30 7.5cm x 2.5cm 10 7.5cm x 5cm	100	1
510	Assorted 4 4cm x 4cm 4 7cm x 2cm 2 7.5cm x 2.5cm	10	1
520	Assorted 8 4cm x 4cm 8 7cm x 2cm 4 7.5cm x 2.5cm	20	1

