

Blue Plaster Range

Whilst incorporating all of the modern features of Dependaplast Washproof Plasters, the Dependaplast Blue Plaster range provides a bright blue colouration that makes them instantly visually detectable should it accidentally come into contact with food.

An aluminium strip located behind the dressing pad is clearly detectable to metal detection devices in food processing facilities.

Dependaplast plasters offer a full and comprehensive range of the most popular shapes and sizes. Plasters are now neatly packed in a variety of boxes, allowing for clear, easy access, and because the boxes are smaller, save on storage too. Ideal for refilling first aid kits or to keep as backup supplies. All Dependaplast plasters are latex-free with a Rayon wound contact layer.



REF	SIZE	BOX QTY	MOQ
558	2cm x 12cm	50	1
448	4cm x 2cm	100	1
541	4cm x 4cm	100	1
449	6cm x 2cm	100	1
544	7.5cm x 2.5cm	100	1
545	7.5cm x 5cm	50	1
543	7cm x 2cm	100	1
549	Anchor	50	1
548	Fingertip	50	1
546	Assorted 30 4cm x 4cm 30 7cm x 2cm 30 7.5cm x 2.5cm 10 7.5cm x 5cm	100	1
542	Assorted 4 4cm x 4cm 4 7cm x 2cm 2 7.5cm x 2.5cm	10	1
540	Assorted 8 4cm x 4cm 8 7cm x 2cm 4 7.5cm x 2.5cm	20	1

