

Clear Plaster Range

Flexible clear PE plastic film is used to ensure a super-thin, flexible, washproof plaster. The membrane is moisture vapour permeable, allowing the skin to breathe without the need for physical perforations. This keeps the dressing securely in place, whilst acting as a complete barrier to water, virus, bacteria, and fungi. High-tech acrylic adhesive provides a really secure fixing, but is skin friendly and low-allergy.

Dependaplast plasters offer a full and comprehensive range of the most popular shapes and sizes. Plasters are now neatly packed in a variety of boxes, allowing for clear, easy access, and because the boxes are smaller, they save on storage too. Ideal for refilling first aid kits or to keep as backup supplies. All Dependaplast plasters are latex-free with a Rayon wound contact layer.



| REF | SIZE | BOX QTY | MOQ |
|------|--|--------------|-----|
| 6930 | 4cm x 2cm | Box of 100 | 10 |
| 6931 | 4cm x 4cm | Box of 100 | 10 |
| 6932 | 6cm x 2cm | Box of 100 | 10 |
| 6933 | 7cm x 2cm | Box of 100 | 10 |
| 6934 | 7.5cm x 2.5cm | Box of 100 | 10 |
| 6942 | 7.5cm x 5cm | Box of 50 | 10 |
| 6939 | Assorted 30 4cm x 4cm 30 7cm x 2cm 30 7.5cm x 2.5cm 10 7.5cm x 5cm | Box of 100 | 10 |
| 6935 | Fingertip | Box of 50 | 10 |
| 6936 | Anchor | Box of 50 | 10 |
| 6940 | Assorted 8 4cm x 4cm 8 7cm x 2cm 4 7.5cm x 2.5cm | Wallet of 20 | 50 |
| 6941 | Assorted 4 4cm x 4cm 4 7cm x 2cm 2 7.5cm x 2.5cm | Wallet of 10 | 50 |
| 6937 | 2cm x 12cm | Box of 50 | 10 |

