

Multisoft Plaster Range

The latest non-woven polyester fabrics are used in Multisoft to provide a soft, highly flexible dressing that is thinner, lighter and more conformable than regular fabric plasters. This product is ideal for activities that require high levels of dexterity, like sewing or PCB assembly. Multisoft uses high-tech acrylic adhesives to provide a really secure fixing, that is skin friendly and low allergy. Multisoft is tough, yet gentle, making it the correct solution for delicate skin, and perfect for children or the elderly.

Dependaplast plasters offer a full and comprehensive range of the most popular shapes and sizes. Plasters are now neatly packed in a variety of boxes, allowing for clear, easy access, and because the boxes are smaller, save on storage too. Ideal for refilling first aid kits or to keep as backup supplies. \ All Dependaplast plasters are latex-free with a Rayon wound contact layer.



REF	SIZE	BOX QTY	MOQ
556	2cm x 12cm	50	1
570	4cm x 2cm	100	1
521	4cm x 4cm	100	1
572	6cm x 2cm	100	1
524	7.5cm x 2.5cm	100	1
525	7.5cm x 5cm	50	1
523	7cm x 2cm	100	1
529	Anchor	50	1
528	Fingertip	50	1
526	Assorted 30 4cm x 4cm 30 7cm x 2cm 30 7.5cm x 2.5cm 10 7.5cm x 5cm	100	1
578	Assorted 4 4cm x 4cm 4 7cm x 2cm 2 7.5cm x 2.5cm	10	1
522	Assorted 8 4cm x 4cm 8 7cm x 2cm 4 7.5cm x 2.5cm	20	1

