

Washproof Plaster Range

Flexible PE plastic film is used to ensure a super-thin, flexible, washproof plaster. The membrane is moisture vapour permeable, allowing the skin to breathe, without the need for physical perforations, keeping the adhesive securely in place, whilst acting as a complete barrier to water, virus, bacteria and fungi. High-tech acrylic adhesive provides a really secure fixing, but is skin friendly and low allergy.

Dependaplast plasters offer a full and comprehensive range of the most popular shapes and sizes. Plasters are now neatly packed in a variety of boxes, allowing for clear, easy access, and because the boxes are smaller, save on storage too. Ideal for refilling first aid kits or to keep as backup supplies. All Dependaplast plasters are latex-free with a Rayon wound contact layer.



REF	SIZE	BOX QTY	MOQ
550	2.2cm	100	1
5311	2.2cm	100	1
557	2cm x 12cm	50	1
530	4cm x 2cm	100	1
531	4cm x 4cm	100	1
532	6cm x 2cm	100	1
534	7.5cm x 2.5cm	100	1
535	7.5cm x 5cm	50	1
533	7cm x 2cm	100	1
539	Anchor	50	1
538	Fingertip	50	1
536	Assorted 30 4cm x 4cm 30 7cm x 2cm 30 7.5cm x 2.5cm 10 7.5cm x 5cm	100	1
554	Assorted 4 4cm x 4cm 4 7cm x 2cm 2 7.5cm x 2.5cm	10	1
553	Assorted 8 4cm x 4cm 8 7cm x 2cm 4 7.5cm x 2.5cm	20	1

