

# UPDATED

# Premier Eye Wash Station

## Now with 500ml Eye Cap Bottles



Robust premier eye wash station.  
Wipe clean, hardwearing, ABS plastic.  
Illuminated, printed instructions.  
Includes integral eye bath.

Dual purpose HSE compliant eye wash station. Wipeable, hard-wearing ABS plastic with illuminated printed instructions; useful if a power cut happens so the station can be located. Dust proof heavy duty cover keeps contents clean. Central Mirror for easy application. Contains easy-to-squeeze bottle with an eye cap, allowing the eye to be irrigated easily and efficiently.



**Eye wash bottle now includes integral eye bath.**

### Premier Eye Wash Station Green/Clear Lid

REF	PRODUCT	SIZE	UNIT	PACK
919	Premier Eye Wash Station	36cmH x 51cmW x 10cmD	1	1
918	Premier Eye Wash Station - empty	36cmH x 51cmW x 10cmD	1	1
933	Refill for Premier Eye Wash Station	-	1	1

#### CONTENTS

2	Dressings, Eye Pad	5	Eye Wash, Pods 20ml
2	Eye Wash, Bottles 500ml	1	Mirror

