

## HYDROFLUORIC ACID

# Treatment Kit for Acid Skin Burns



Even in a diluted form, hydrofluoric acid is highly corrosive and toxic. It can severely damage the skin and eyes, causing severe burns which are extremely painful and can leave long-term scarring with bone and tissue damage.

The dedicated Hydrofluoric Acid Treatment Kit contains all the critical items required to treat Hydrofluoric Acid burns, such as Calcium Gluconate Gel which helps to neutralise the powerful fluoride ion, and sterile wound dressings to reduce the risk of infection; whilst providing a suitable level of protection for the first aid responder.



*The Hydrofluoric Acid Treatment Kit includes chemical resistant nitrile gloves to ensure the safety of the First Aid Responder.*



Tested for Microorganism Hazards



Tested for Watertightness and Low Chemical Protection



### Hydrofluoric Acid Treatment kit

#### CONTENTS

2	Bandages, Conforming 7.5cm x 4m	2	Gloves, Nitrile Pairs
2	C-Gel Calcium Gluconate, 40g	1	Shears
2	Dressings, Low Adherent 10cm x 10cm		

5975 Hydrofluoric Acid Treatment Kit for Acid Skin Burns  
in Red/Orange Compact Aura  
20.5cmH x 20cmW x 6.5cmD

