

'High quality  
CPR and rapid  
defibrillation  
save lives'



Fast and easy set-up buttons for language, scenarios, volume, and play/pause.



Exclusive to the Prestan AED Trainer is the Pad Sensing System, an embedded 'switch' in the pad itself. The switch 'senses' that the AED Trainer pads are attached to the manikin torso.



Optional Trainer Remote Control.

PRESTAN®

## AED Trainer

The feature-packed Prestan Professional AED Trainer offers a number of attributes usually found on more expensive models.

The Prestan Professional AED Trainer provides quality, durability, and affordability. Clear, confident voice instructions provide easy-to-follow directions, plus realistically configured scenarios that makes training students easier.

View the full range of Prestan manikins and video tutorials now at [www.prestan.co.uk](http://www.prestan.co.uk)



# PRESTAN® AED Trainer

- Prestan quality, durability and affordability.
- Unique pad sensing system.
- Revolutionary software and language module.
- Instructor friendly setup and use.
- Languages available: English/French and English/Spanish (more to follow).
- Exclusive to the Prestan AED Trainer is the Pad Sensing System, an embedded 'switch' in the pad itself.
- The switch 'senses' that the AED Trainer pads are attached to the manikin torso.
- Pads are already connected and may be used at least 25-30 times during CPR/AED classes.
- The Pads are made of the same foam material used on functioning AEDs and are compatible with all CPR manikins on the market.
- The silicone based adhesive on the pads, coupled with the foam structure, ensures easy application and removal.
- Clear and calm voice prompts.
- Chest compression metronome (recording of the actual clicker in a Prestan Professional Adult Manikin).
- Ascending counts in 10's, ending at 30.
- Made of 'super-strong', durable polycarbonate.
- Friction hinges allow for smooth closing of the lid.
- Integrated magnets hold the lid securely in place.



Contact [sales@reliancemedical.co.uk](mailto:sales@reliancemedical.co.uk) or 08456 448808 and we will tell you how to contact your local distributor.

The Prestan Professional AED Trainer provides excellence, sturdiness, and economy. Distinct voice commands deliver assertive, easy-to-follow directions, including predetermined scenarios that make teaching students stress-free and much easier. The unit includes one set of adult training pads, batteries, instructions for use and a carrying case. Remote control is available as a separate option. One Year Warranty.



PRESTAN AED TRAINER

7440	Prestan AED Trainer English/Spanish Single	Unit 1	MOQ 1
7442	Prestan AED Trainer English/French Single	Unit 1	MOQ 1



PRESTAN AED TRAINER

7441	Prestan AED Trainer English/Spanish Pack of 4	Unit 4	MOQ 1
7443	Prestan AED Trainer English/French Pack of 4	Unit 4	MOQ 1



PRESTAN AED TRAINER REMOTE CONTROL

7444	Prestan AED Trainer Remote with 2 AA Batteries Single	Unit Single	MOQ 1
7445	Prestan AED Trainer Remote with 2 AA Batteries Pack of 4	Unit 4	MOQ 1



PRESTAN AED TRAINER LANGUAGE MODULE

7446	Prestan AED Trainer Language Module English/Spanish Single	Unit Single	MOQ 1
7447	Prestan AED Trainer Language Module English/Spanish Pack of 4	Unit 4	MOQ 1
7448	Prestan AED Trainer Language Module English/French Single	Unit Single	MOQ 1
7449	Prestan AED Trainer Language Module English/French Pack of 4	Unit 4	MOQ 1



PRESTAN AED TRAINER PADS

7450	Prestan AED Trainer Pads Single Pair	Unit Single	MOQ 1
7451	Prestan AED Trainer Pads Pack of 4 Pairs	Unit 4	MOQ 1

