

'High quality CPR and rapid defibrillation save lives'



*Fast and easy set-up buttons for language, scenarios, volume, play/pause, and child option.*



*Exclusive to the Prestan AED Trainer is the Pad Sensing System; an embedded 'switch' in the pad itself. The switch 'senses' that the AED Trainer pads are attached to the manikin torso.*

View the full range of Prestan manikins and video tutorials at [www.prestan.co.uk](http://www.prestan.co.uk)

## PRESTAN® AED UltraTrainer™

CPR feedback and coaching.

Child defibrillation.

Pad Sensing System.

3-year warranty.

The Prestan AED UltraTrainer™ offers features not found on any other universal AED trainer, and is the most accurate representation of how current AEDs in the marketplace perform. This compact and lightweight unit will give students a realistic training experience with training pads that sense when they are placed on a manikin. In addition, the pads are pre-connected and voice prompts are clear and calm.

Instructors can customise their training classes with options to turn the compression metronome on or off, and/or the "Give Breaths" voice prompt. With a 3-year warranty and durable long-lasting pads, the AED UltraTrainer™ is built to last.



# AED UltraTrainer™

The Prestan AED UltraTrainer™ is a compact, lightweight solution with full customisation options. Each unit includes one set of adult/child training pads, batteries, instructions for use, a carrying case, and 3-year warranty.

**Prestan is the First and only Universal AED Trainer with:**

- Dual Graphic adult / child, training pads.
- Pad Sensing System.
- Pre-connected pads.
- Child training button.
- Child AED specific voice prompts.
- Automatic or Semi-automatic shock simulation.
- CPR coaching, "Give Breaths" voice prompt.
- Turn compression metronome on or off.

- Turn on or off "Give Breaths" prompts.
- Low battery warning.
- 3-year warranty.

**Other features:**

- Two languages.
- 2015 Guidelines compliant and upgradeable.
- Play / pause button.
- Five training scenarios.



**Prestan AED UltraTrainer English/Spanish**

REF	PRODUCT	SIZE	UNIT	PACK
7660	Prestan AED UltraTrainer English/Spanish	18.5cmH x 13cmW x 5cmD	1	1

CONTENTS				
1	AED UltraTrainer™	1	Instructions	
1	Adult/Child Pads	1	Connector Buttons	
1	Batteries	1	Carry Case	

**Prestan AED UltraTrainer English/Spanish**

REF	PRODUCT	SIZE	UNIT	PACK
7661	Prestan AED UltraTrainer English/Spanish	18.5cmH x 13cmW x 5cmD	Pack of 4	1

CONTENTS				
4	AED UltraTrainer™	4	Instructions	
4	Adult/Child Pads	4	Connector Buttons	
4	Batteries	4	Carry Cases	



Contact [sales@reliancemedical.co.uk](mailto:sales@reliancemedical.co.uk) or 08456 448808 and we will tell you how to contact your local distributor.

View the full range of Prestan manikins and video tutorials at [www.prestan.co.uk](http://www.prestan.co.uk)



**Prestan AED UltraTrainer Adult/Child Replacement Training Pads**

REF	PRODUCT	SIZE	UNIT	PACK
7662	Prestan AED UltraTrainer Adult/Child Replacement Training Pads	Adult/Child	Set of 2	1
7663	Prestan AED UltraTrainer Adult/Child Replacement Training Pads	Adult/Child	Set of 8	1

