

Bandages

Religrip Elasticated Tubular Bandages



Supportive, comfortable, elasticated tubular bandage. Ideal treatment for soft tissue injuries, sporting injuries, dislocations and sprains. Religrip provides an even radial support and pressure for weakened joints and muscles, also aiding rehabilitation after the removal of plaster of paris casts. It remains in place without ties or tapes.

Size Chart

| Size | Width | Application |
|--------|--------|---|
| Size A | 4.5cm | Child's small limb |
| Size B | 6.25cm | Child's medium limb Adult small arm |
| Size C | 6.75cm | Child's large limb Adult medium arm Adult small ankle |
| Size D | 7.5cm | Adult large arm Adult medium ankle Adult small leg |
| Size E | 8.75cm | Adult large ankle Adult medium leg Adult small thigh |
| Size F | 10cm | Adult large leg Adult medium thigh |
| Size G | 12cm | Adult large thigh |

NB: The chart above is a guidance only.

Religrip Elasticated Tubular Bandage

Natural

Supportive, comfortable, radial bandage. Ideal for soft tissue injuries, dislocations and sprains. Natural colour for general use. Range of sizes for all parts of the body. Available in seven widths.

| Code | Product Name | Size | Unit | Pack |
|------|--------------------------------------|------------|---------|----------|
| 461 | Religrip Elasticated Tubular Bandage | Size A 10m | Unit: 1 | Pack: 1 |
| 462 | Religrip Elasticated Tubular Bandage | Size B 10m | Unit: 1 | Pack: 1 |
| 463 | Religrip Elasticated Tubular Bandage | Size C 10m | Unit: 1 | Pack: 1 |
| 464 | Religrip Elasticated Tubular Bandage | Size D 10m | Unit: 1 | Pack: 1 |
| 465 | Religrip Elasticated Tubular Bandage | Size E 10m | Unit: 1 | Pack: 1 |
| 466 | Religrip Elasticated Tubular Bandage | Size F 10m | Unit: 1 | Pack: 1 |
| 467 | Religrip Elasticated Tubular Bandage | Size G 10m | Unit: 1 | Pack: 1 |
| 471 | Religrip Elasticated Tubular Bandage | Size A 1m | Unit: 1 | Pack: 10 |
| 472 | Religrip Elasticated Tubular Bandage | Size B 1m | Unit: 1 | Pack: 10 |
| 473 | Religrip Elasticated Tubular Bandage | Size C 1m | Unit: 1 | Pack: 10 |
| 474 | Religrip Elasticated Tubular Bandage | Size D 1m | Unit: 1 | Pack: 10 |
| 475 | Religrip Elasticated Tubular Bandage | Size E 1m | Unit: 1 | Pack: 10 |
| 476 | Religrip Elasticated Tubular Bandage | Size F 1m | Unit: 1 | Pack: 10 |
| 477 | Religrip Elasticated Tubular Bandage | Size G 1m | Unit: 1 | Pack: 10 |

