

Relipad

Low Adherent Dressing Pads



Sterile low-adherence absorbent pads.

Double-sided; prevents misapplication.

For minor abrasions and lacerations.

Sterile, low adherent dressing pads, designed to be applied to a wound from either side, thus preventing incorrect application. The dressing is low adherent with high absorbency, effectively drawing exudate away from the wound. For use on abrasions, lacerations and moderately exuding wounds. Soft, comfortable and conformable. Secure with tape or a conforming bandage.



REF	PRODUCT	SIZE	UNIT	PACK
381	Relipad Low-Adherent Dressing Pads Sterile	5cm x 5cm	Box of 100	1
382	Relipad Low-Adherent Dressing Pads Sterile	7.5cm x 7.5cm	Box of 100	1
383	Relipad Low-Adherent Dressing Pads Sterile	10cm x 10cm	Box of 100	1

