

TraumaFix® MILITARY

Designed for major bleed trauma



Fabric waterproof packaging.

Pressure bandage.

Easy and fast self administration.

Applies pressure whilst the casualty is being moved.

Green camouflage packaging, with skin camouflage bandage.

What do you have in your First Aid Kit to deal with a major bleed injury?

Ambulance dressings have their place, but they are not substantial enough to cope with a major trauma bleed, and applying enough pressure with a light conforming bandage will be near impossible. Abdo-pads might be substantial enough to cope with the bleed, but they need a suitable pressure bandage to attach them, which then needs to be secured. This is very time consuming, difficult with one injury, dangerous with multiple bleeds.

All too often, non-sterile sanitary towels are used as major bleed dressing pads. Typically these contain super-absorber powders that clot blood into a gel on contact. Just one grain of this powder has the capacity to cause a deadly blood clot if it were to ingress into the casualties venous system.

TraumaFix® offers a substantial, highly absorbent dressing pad, capable of holding over ten times its own weight in fluid. A strong bandage allows direct pressure to be applied quickly and maintained easily. This is secured instantly with a Velcro® style fastening.

Multiple major bleed injuries can be stemmed quickly, allowing the casualty to be moved easily. The dressing is a significant advance in the treatment of trauma injuries.

REF	PRODUCT	DESCRIPTION	SIZE	UNIT	MOQ
1964	TraumaFix Military Dressing	Sterile	10cm x 19cm	1	10
1965	TraumaFix Military Dressing	Sterile	20cm x 19cm	1	10

**IMPROVED
NOW WITH
PRESSURE
PILLOW**



Innovative new pillow increases the pressure to stem major blood flow.

Major bleed technology.

Amazing fluid absorption.



Place highly absorbent dressing pad with pressure pillow over trauma...



Direct pressure is applied with pressure pillow and 4.5m of pressure bandage...



The bandage is secured instantly with a quick velcro style fastening...



Major bleeding is stemmed quickly, allowing the casualty to be moved easily.

