

Will your AED be there when you need it?

Secure and Protect your lifesaving AED's with our alarmed, AED Storage Cabinet



- Toughened durable glass door.
- Compatible with all AED models.
- Safe and visible.
- Alarm sounds when opened.
- Wall mountable.

The AED wall cabinet sounds an alarm instantly once the door has been opened. The alarm alerts people in the vicinity of the emergency, improves security and acts as a deterrent against unnecessary tampering. Compatible with all major AED models.

The cabinet is constructed from a white powder coated metal, with a stone enameled finish and toughened glass door for durability and visibility.

Secure your vital AED's, provide rapid access, maintain in dry, dust-free conditions, and be ready for an emergency with the AED wall mounted alarmed cabinet.



Storage Cabinet

REF	PRODUCT	SIZE	UNIT	PACK
2102	AED Alarmed Metal Storage Cabinet White/Glass/Alarmed	Internal: 43.4cmH x 43.4cmW x 21.7cmD External: 44.4cmH x 44.4cmW x 22.7cmD	1	1

