

Heated Outdoor Metal Wall Cabinet*

to secure and protect your lifesaving AED's

- Compatible with all AED's
- Multiple Security Options
- AED Signage
- Thermostat Controlled
- Internal LED Light

- **Mechanical Keypad**

The mechanical keypad has a permanent blue LED light which allows the user to access their AED with their chosen numeric combination.

- **Digital Temperature and Thermostat Display**

The digital temperature display has an illuminated LCD screen displaying the internal temperature. The user can view the display, set the temperature and be certain that the defibrillator is being stored at the correct temperature.

- **Integrated Alarm with Remote Control**

Once activated, the integrated alarm will sound when the cabinet is opened. The remote control allows you to deactivate the alarm prior to opening the cabinet.

- **Manufacturers 12 Month Warranty**

* AED not included



Our AED wall mounting cabinet offers you security, practicality, high visibility, outdoor weather resistance and peace of mind in the event of an emergency.

The cabinet is manufactured from thick carbon steel with the addition of a shatterproof transparent window, allowing you to quickly check the status of your defibrillator. Finished in visible yellow, with the added features of a waterproof seal, user regulated temperature settings and a magnetic LED light which is operated when the door opens. The cabinet also includes AED instruction icons to guide users on the correct procedure.



COMPONENTS

Mechanical keypad
Alarm (120DB)
Heating system
Thermostat
Magnetic LED light
9V blue LED light
9V transformer

SPECIFICATION

Weight: 11.5kg
Size: 50cm x 40cm x 22cm
Internal size:
37cm x 38cm x 17cm
Material: 1.5mm Steel
IP56 Waterproof rating
Option to hard wire

POWER SUPPLY

Main power: 220V
Alarm: 9V Battery
Thermostat/Heating system:
220V, 100W, UK 3 pin plug
Magnetic LED light: 220V, 5W
Alarm: 9V battery
Blue LED illumination:
9V, 0.5W

Outdoor Metal Wall Cabinet with Integrated Alarm

REF	PRODUCT	SIZE	UNIT	PACK
2105	AED Alarmed, Outdoor, Heated, Metal Cabinet Yellow/Alarmed/Heated	Internal: 50cmH x 40cmW x 22cmD External: 37cmH x 38cmW x 17cmD	1	1